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3. FMX-energiemodule
4. FMX chirurgische instrumenten
5. FMX-voetschakelaar

English

Introduction	3	Generator Start-Up Instructions.....	5
Indications for Use.....	3	Foot Pedal Installation.....	5
Contraindications	3	Power Module and Surgical Instrument Installation	5
Important Information	3	Preparing Generator for Use	6
Cautions	3	Operating Instructions.....	6
Precautions.....	3	Cleaning Instructions.....	7
Warnings	4	Generator Power-Down Instructions	7
Front Panel Layout.....	4	Symbols.....	7
Rear Panel Layout.....	4	Specifications	7
Generator Installation.....	5	Troubleshooting	9
Generator Preparation.....	5	Warranty.....	10

Introduction

Domain Surgical's FMX Ferromagnetic Surgical System is an intelligent surgical platform that provides precise control of surgical thermal energy, without passing any electrical current through the patient. The FMX system consists of a generator, power module, and single-use surgical instruments.

Indications for Use

The FMX Ferromagnetic Surgical System is indicated for cutting and coagulation of soft tissue. For specific indications of individual FMX surgical instruments, see their respective Directions for Use.

Contraindications

For specific contraindications of individual FMX surgical instruments, see their respective Directions for Use.

Important Information

Prior to use of the FMX system, ensure that all package inserts, including warnings, precautions, and instructions are read and understood.

Cautions

- Federal (US) law restricts this device to sale by or on the order of a physician.
- This manual and the equipment described within are for use only by qualified medical personnel trained in the surgical procedures to be performed.
- Prior to initial use of the FMX system, ensure that all package inserts, including warnings, precautions and instructions are read and understood.


Precautions

- The FMX system and components should only be used by qualified medical personnel trained in the surgical procedures to be performed.
- Take care when handling the surgical instruments to prevent possible damage to the tip.
- Do not use abrasive materials or metal instruments to clean the surgical tip. This can damage the tip and cause the device to fail to function as intended. Gauze moistened with water or saline is recommended for removing heavy coagulum.
- Use the lowest power level and the shortest activation time necessary to achieve the desired surgical effect.
- Prior to use, carefully inspect the FMX surgical instruments for any defects. DO NOT USE if the packaging, cable insulation or the connector appear damaged.


Warnings

- ⚠ Read all instructions before use.
- ⚠ Use ONLY FMX surgical instruments with the FMX Generator.
- ⚠ DO NOT reuse, re-sterilize or reprocess "Single Use Only" labeled accessories. Reuse of "Single Use Only" devices risks cross-contamination of substances, viable bacteria, and viruses between patients.
- ⚠ DO NOT use in the presence of flammable liquids, gasses or materials. DO NOT activate the handpiece until vapors from alcohol-based skin prep agents have dissipated.
- ⚠ Flammable agents used for cleaning or disinfecting, or as solvents of adhesives, should be allowed to evaporate before the application of high frequency surgery. There is a risk of pooling of flammable solutions under the patient or in body depressions such as the umbilicus, and in body cavities such as the vagina. Any fluid pooled in these areas should be mopped up before high frequency surgical equipment is used.
- ⚠ Attention should be called to the danger of ignition of endogenous gases. Some materials, for example cotton and gauze, when saturated with oxygen may be ignited by heat produced in normal use of the high frequency surgical equipment.
- ⚠ Be sure to allow sufficient time for surgical tip to return to ambient temperature to avoid unintended burns to patient tissue or the surgical staff.
- ⚠ DO NOT place the instruments on the patient drape. The tip can remain hot after deactivation and can potentially melt through the drape. Placement of the instrument on a wet, sterile towel when not in use will prevent the possibility of melting the drape if the tip remains hot or if the instrument is inadvertently activated with the foot pedal.
- ⚠ DO NOT inadvertently press the foot pedal. In order to prevent the instrument from being unintentionally activated, place the foot pedal in a location necessitating deliberate action in order to activate the generator.
- ⚠ DO NOT bend the tip of the surgical instruments. Care should be taken not to deliberately bend the tip as damage will prevent proper function of the handpiece.
- ⚠ Only the surgeon using the instrument should operate the foot pedal. Unintentional activation may occur if the foot pedal is activated by a separate user, which may result in patient or user injury.
- ⚠ Smoke generated during surgical procedures can be irritating and potentially harmful to surgical personnel. It is recommended to use surgical masks and adequate ventilation or evacuation of smoke.
- ⚠ ELECTRICAL SHOCK HAZARD: Do not remove the cover of the FMX Generator. Refer to authorized personnel for service.
- ⚠ To allow for proper cooling, the generator should not be placed in an enclosed area while in operation. Appropriate cooling is necessary to avoid overheating of the generator.
- ⚠ The power cord supplied should be connected to a properly grounded receptacle. Extension cords and/or adapter plugs should not be used. If necessary, the potential equalization terminal may be used to establish electrical equipotentiality with the grounding bus bar of the electrical installation or with other equipment (per the requirements of IEC 60601-1).
- ⚠ Interference produced by the operation of the FMX Generator may adversely influence the operation of other electronic equipment.
- ⚠ Operating room staff should never contact the surgical tip while the generator is activated or while the tip is cooling down, as injury may result.
- ⚠ When high frequency surgical equipment and physiological monitoring equipment are used simultaneously on the same patient, any monitoring electrodes should be placed as far as possible from the surgical electrodes. Needle monitoring electrodes are not recommended. In all cases, monitoring systems incorporating high frequency current limiting devices are recommended.
- ⚠ Avoid contact with implanted electrical lead wires. The protective thermoplastic insulation on some wires may melt at high temperatures.
- ⚠ DO NOT clean the tip of the surgical instrument with an abrasive pad or metal instrument.
- ⚠ Failure of the FMX Generator could result in an unintended increase of output power.
- ⚠ Ensure that the sound volume on the generator is adequately adjusted so that the activation tones can be clearly heard. The activation tone is intended to alert the user that the device is active. This will help to prevent unintended activation of the device, which could result in patient or user injury.
- ⚠ Ensure that the LCD touchscreen display and audio speaker are functioning when turning on the generator. Several audible beeps should be heard during the start-up process, and a splash screen will be displayed.
- ⚠ Ensure that all power cords and connector cables are properly organized and placed so as to not create a tripping hazard.
- ⚠ No modification of this equipment is allowed.
- ⚠ Covers or access panels are to be removed only by authorized service personnel.

Front Panel Layout

Rear Panel Layout

Generator Installation

Unpack and inspect the generator for any signs of physical damage to the front panel or housing. Make sure that there are no exposed wires or other damage that would allow the internal electrical components to be accessed.


The generator may be placed on a cart or on any sturdy table, shelf, or platform.

NOTE: Adequate space surrounding the generator is essential to allow air flow for proper cooling of the generator. Refer to hospital procedures or local codes for additional requirements.


NOTE: The power cord supplied should be connected to a properly grounded receptacle. Extension cords and/or adapter plugs should not be used.

Generator Preparation


Plug the provided power cord in to the mains power cord connection on the rear panel of the generator.

Plug the power cord into a grounded receptacle. Activate the rear panel power button on the back panel of the generator.

When the generator is plugged in and the rear panel power button on the back panel is switched 'on', the front panel power button will blink blue. The LCD touchscreen display on the front of the generator will not illuminate until the front panel power button is pressed.


Generator Start-Up Instructions

Pressing the **front panel power button** will activate the unit and start a series of internal self-tests. These tests will take a few seconds and are designed to assure that the generator is ready to operate as intended.

If any of the self-tests fail, a message will be displayed on the touchscreen. These error messages indicate that the generator requires attention. See the "Troubleshooting" section of this document for the proper course of action if this occurs.


When the generator starts without a power module and surgical instrument connected, a screen indicating that no device is connected will be displayed, and the speaker will be set to the minimum volume.

Any alerts or alarms during operation will be displayed on the touchscreen and accompanied by several audible beeps.

Foot Pedal Installation

If desired, plug the foot pedal into the foot pedal connection on the back panel of the generator.


CAUTION: Use only the FMX Foot Pedal with the FMX Generator.

Power Module and Surgical Instrument Installation


The FMX System requires the use of an FMX Power Module to manage the delivery of energy to the surgical instruments.

FMX power modules are available in either resterilizable or single-use disposable configurations. The FMX Reusable Power Module has a dark blue cable and should have been sterilized by your facility. The FMX Disposable Power Module has a white cable and is packaged in a sealed sterile pouch. Both power modules are compatible with all FMX surgical instruments.

Connect the power module to the surgical instrument by inserting it into the instrument and push until you feel a "click".


NOTE: If using an FMX Reusable Power Module, an air cooling line (which was packaged with individual FMX surgical instruments) must be attached to both ends of the power module cable by firmly inserting and twisting the tubing into the black receptacle.


NOTE: The FMX Disposable Power Module features an integrated air cooling line, and does not require connection of the additional air line.

Connect the power module to the generator by inserting it into the connection port on the front of the generator.

An audio indicator will sound, and the LCD Touchscreen Display will display 2 green checkmarks, indicating that both the power module and surgical instrument have been properly connected to the generator.

NOTE: If using an FMX Reusable Power Module, the generator screen will display how many procedures remain in the service life of the power module immediately after the power module is connected. If using an FMX Disposable Power Module, the generator screen will display "Single Use Only."

NOTE: For detailed installation instructions for specific FMX surgical instruments, please see their respective Directions for Use.


CAUTION: The technology employed in the FMX System is unique. Electrosurgical handpieces from other manufacturers are not compatible and may introduce unintended hazards or damage to the generator if they are used with the FMX Generator.

Preparing Generator for Use

The default display is the power control display. This screen displays the current power settings, an icon indicating current audio level setting, an icon indicating if a foot pedal is connected, and an icon indicating which instrument (if any) is connected.

Some FMX instruments feature dual activation buttons, which allow the surgeon to quickly toggle between user-selected high and low power settings. When an FMX instrument with dual activation buttons is connected to the generator, the power control display will display 2 separate power settings. Pressing the up or down arrow buttons

on the generator display will change the power level until the arrow is released or the maximum/minimum power level is reached.

NOTE: Set the power level to the lowest power level that provides the intended surgical effect.

NOTE: The power control screen varies slightly depending on what surgical instrument is connected. Refer to the each instrument's Directions for Use for specific information.

Touching the 'speaker' icon in the upper left of the power control display brings up the audio control display.

The default volume setting is the lowest possible volume level. Pressing the up and down arrow buttons on the audio control display will change the volume of the tones that are sounded in normal operation. Touching the 'back' arrow will return the display to the power control display. If there is no activity on the audio control display for 5 seconds, the display will automatically return to the power control display.

Operating Instructions

There are two ways to activate the surgical instrument.

- Completely press and hold either of the activation buttons on the surgical instrument to activate the selected power setting.
- Step on the desired foot pedal to activate the selected power setting. The dark blue pedal is used to activate the high power setting, and the light blue pedal is used to activate the lower power setting.

To deactivate the instrument, both the button on the instrument and the foot pedal (if attached) must be released.

CAUTION: Be sure to allow sufficient time for heating element on the surgical tip to return to ambient temperature before using as a probe/retractor. The activation tone is audible for a few seconds following release of the activation button to indicate that the surgical tip may still be hot.

CAUTION: Explosion and Fire Hazard: Because the FMX instruments can produce temperatures high enough to ignite some flammable gasses, do not use in the presence of flammable anesthetic mixtures with air, oxygen or nitrous oxide.

Cleaning Instructions

The generator, mains power cord and foot pedal are cleanable but are not intended to be sterile. Do not allow fluids to enter the housing of the generator or foot pedal.

- Before cleaning, unplug any power module that may be attached to the generator, and unplug the generator from the power source.
- Do not immerse the generator or foot pedal in any liquid.
- Do not use an abrasive cloth or abrasive cleaner anywhere on the generator.

The generator, power cord and foot pedal may be cleaned with a soft cloth dampened with water, alcohol, mild soaps, detergents or disinfectants. Non-flammable agents should be used for cleaning and disinfection wherever possible.




Generator Power-Down Instructions

To power-down the generator, press the front panel power button and hold for 1 second. A short beep will sound, and the message "Powering Down" will display on the LCD touchscreen.

When power-down is complete, the screen will go dark. It is now safe to turn off the rear power button and disconnect the power cord.

Symbols

The following symbols are found on the front and back panels of the generator as well as on the accompanying documents of the generator and accessories.

	Consult Instructions For Use
	Single-use Only
	DO NOT USE if package appears damaged
	Serial Number
	Lot Number
	Date of Manufacture
	Type CF Applied Part
	Warning: Dangerous Voltage
	Non-ionizing Radiation
	Temperature Limitation
	Caution
	Alternating Current
	Conductor may be used to provide potential equalization between the equipment and the installation busbar

Specifications

Operating Environment Conditions

- Temp: 10° to 30° C (50° to 86° F)
- Humidity: 5% to 85% Non-condensing

Storage Conditions

- Temp: -20° to 65° C (-4° to 149° F)
- Humidity: 15% to 85% Non-condensing

Shipping Conditions

- Temp: -20° to 65° C (-4° to 149° F)
- Humidity: 15% to 85% Non-condensing

Moisture Protection

- Generator: IPX1
- Foot Pedal: IPX7 minimum

Generator Dimensions

- 14" x 4" x 12"
(35.6 cm x 10.2 cm x 30.5 cm)

Generator Weight

- Approximately 20 lbs. (9 kg)

Power Requirements

- Voltage: 100/110/220/230 Vac
- Current: 4/2A rms max
- Frequency: 50/60 Hz

Power Output

- 60W Max

Fuses

- For 100/110V: T 5A L 250V
- For 220V/230V: T 2.5A L 250V

Fundamental Output Frequency

- 40.68 MHz

Mode of Operation

- Intermittent Operation
10 sec. ON, 30 sec. OFF

Protection Against Electric Shock

- Class I

Applied Part Classification

- Type CF

Guidance and manufacturer's declaration – electromagnetic emissions

The FMX Ferromagnetic Surgical System is intended for use in the electromagnetic environment specified below. The customer or the user of the FMX Ferromagnetic Surgical System should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The FMX Ferromagnetic Surgical System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The FMX Ferromagnetic Surgical System is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded:
Harmonic emissions IEC 61000-3-2	Class A	WARNING: This equipment/system is intended for use by healthcare professionals only. This equipment/ system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting or relocating the FMX Ferromagnetic Surgical System or shielding the location.
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration

The FMX Ferromagnetic Surgical System is intended for use in the electromagnetic environment specified below. The customer or the user of the FMX Ferromagnetic Surgical System should assure that it is used in such an environment.

Immunity test	IEC 60601 test Level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV line to line ±2kV line to earth	±1kV line to line ±2kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	>95% dip in mains voltage value for 0.5 cycle 60% dip in mains voltage value for 5 cycles 30% dip in mains voltage value for 25 cycles >95% dip in mains voltage value for 5 seconds	100% reduction for 0.5 cycle 60% reduction for 5 cycles 30% reduction for 25 cycles 100% reduction for 250 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the FMX Ferromagnetic Surgical System required continued operation during power mains interruptions, it is recommended that the system be powered from an uninterruptable power supply or battery.
Power frequency (50/60 Hz) magnetic Field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Recommended separation distances between portable and mobile RF communications equipment and the FMX Ferromagnetic Surgical System

The FMX Ferromagnetic Surgical System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the system as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d=[1.2]√P	80 MHz to 800 MHz d=[1.2]√P	800 MHz to 2.5 GHz d=[2.3]√P
0.01	.12	.12	.23
0.1	.38	.38	.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – electromagnetic immunity

The FMX Ferromagnetic Surgical System is intended for use in the electromagnetic environment specified below. The customer or the user of the system should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	3Vrms	Portable and mobile RF equipment should be used no closer to any part of the FMX Ferromagnetic Surgical System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency transmitter. Recommended separation distance: d=[1.2]√P 80 MHz to 800 MHz: d=[1.2]√P 800 MHz to 2.5 GHz: d=[2.3]√P Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3V/m 80 MHz to 2.5 GHz	3V/m	

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Troubleshooting

Alert Message	Action Required
HANDPIECE NOT FOUND	Verify that the surgical instrument is securely connected to the generator.
HANDPIECE BROKEN	Disconnect and reconnect handpiece to the generator. If problem persists, replace with new handpiece.
HANDPIECE COM FAIL	
TIP/HANDPIECE FAIL	
TIP NOT CONNECTED	Verify that the surgical tip is securely connected to the handpiece.
TIP ERROR	Disconnect the surgical tip from the handpiece and reconnect the tip to the handpiece. If problem persists, replace with new tip.
TIP NOT INITIALIZED	Disconnect the surgical tip from the handpiece and replace with a new tip. If problem persists, contact Domain Surgical for assistance.
POWER MODULE BROKEN	Disconnect and reconnect power module to the generator. If problem persists, replace with new power module.
TOOL/PM FAIL	
POWER MODULE COM FAIL	
TOOL NOT CONNECTED	Verify that the surgical instrument is securely connected to the power module.
TOOL ERROR	Disconnect the surgical instrument from the power module and reconnect the instrument to the power module. If problem persists, replace with new instrument.
TOOL NOT CONFIGURED	Disconnect the surgical instrument from the power module and replace with a new instrument. If problem persists, contact Domain Surgical for assistance.
RELEASE TOOL BUTTON	Release the activation button on the surgical instrument. If problem persists, disconnect the surgical instrument from the power module and reconnect the instrument to the power module. If problem persists, replace with new instrument.
'n' RPM USES REMAINING	This warning indicates that the reusable power module has less than 5 uses remaining in its service life.
PM USE COUNT EXCEEDED	The service life of the power module has been exceeded. Replace with new power module.
AIR LINE BLOCKED	Uncoil the handpiece cable to ensure there are no kinks. If problem persists, replace with new handpiece.
AIR PUMP FAIL	Disconnect the handpiece from the generator. Uncoil the handpiece cable to ensure there are no kinks. Reconnect handpiece to generator. If problem persists, contact Domain Surgical to return the generator for service.
RPM AIR FLOW FAILURE	Verify that an air line is securely attached to both ends of the power module. If problem persists, contact Domain Surgical for assistance.
GENERATOR TOO HOT	Leave the generator plugged in. Power-down the generator using the Front Panel Power Button only. Allow the generator to cool until the internal cooling fans turn off. If the error persists, contact Domain Surgical to return the generator for service.

Alert Message	Action Required
PRE-AMP I TOO HIGH	Release the activation switch on the handpiece and/or foot pedal to clear the alarm. If the error persists, contact Domain Surgical to return the generator for service.
RF-AMP I TOO HIGH	
RF BIAS ERROR	
RF OUTPUT TOO HIGH	
POWER CONTROL ERROR	
FOOT SWITCH ERROR	Make sure that the foot pedal is not depressed. Unplug and plug in the foot pedal to the rear of the generator. If the error persists, contact Domain Surgical to return the generator and foot pedal for service.
NO SWR METER	Release the activation switch on the handpiece and/or foot pedal to clear the alarm. Completely power-down the generator (both Front Panel Power Button and Rear Panel Power Button). Power-up the generator. If the error persists, contact Domain Surgical to return the generator for service.
RF 50 OHM FAIL	Completely power-down the generator (Rear Panel Power Button). Power-up the generator. If the error persists, contact Domain Surgical to return the generator for service.
RF 100 OHM FAIL	
RF OPEN FAIL	
SOFTWARE ERROR	
CRC FAIL - PROGRAM	
CRC FAIL - LOG	
DATA FAILURE	
POWER ON FAIL	
PS 24V FAIL	
TEMP SENSOR FAIL	
PRESS. SENSOR FAIL	
LCD HEARTBEAT LOST	
CPU FAILURE	
AC FAIL CHECK POWER CABLE	Completely power-down the generator (Rear Panel Power Button). Verify that the power cable is connected properly. Power-up the generator. If the error persists, contact Domain Surgical to return the generator for service.
CALIBRATION NEEDED	Contact Domain Surgical to return the generator for service.
CRC FAIL - CONFIG	
CRC FAIL - CALIB	Contact Domain Surgical for assistance.
UNSUPPORTED TIP	
UNSUPPORTED ACCESSORY	
UNSUPPORTED TOOL	

Please be sure to retain any faulty surgical instruments and return to Domain Surgical for evaluation. Before returning a generator, power module, or instrument, please call Domain Surgical at (801) 924-4959 to obtain a Return Material Authorization number and shipping instructions.

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Notes

For customer support inquiries,
please contact us at 888.666.4484

Domain Surgical, Inc.
1370 South 2100 East
Salt Lake City, UT 84108
888.666.4484
omni-guide.com

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Domain Surgical, FMX, FMwand, and FMsealer are registered trademarks of Domain Surgical, Inc. All content is copyright of Domain Surgical, Inc. © 2017. Visit www.omni-guide.com/patents for information about the coverage of these products by U.S. and/or foreign patent rights. LB-001 Rev. 0

☎ 888.666.4484 | 617.551.8444
✉ customerservice@omni-guide.com
🌐 www.omni-guide.com