**OMNIGUIDE ADVANCED ENERGY FIBERS:
BEAMPATH OTO-M
INSTRUCTIONS FOR USE**

OmniGuide, Inc.
4 Maguire Rd.
Lexington MA 02421 USA
Phone: (617) 551-8444 Fax: (617) 551-8445
www.Omni-Guide.com

Caution: Federal Law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

SECTION I: DEVICE DESCRIPTION

(Humidity limitation: The range of humidity to which the fiber can be safely exposed)

The OmniGuide CO₂ Laser Fibers enable the delivery of CO₂ laser energy from a laser source to the target tissue. The fiber attaches to an OmniGuide Fiber Enabled Laser System. Along with the CO₂ laser energy, the BeamPath-M series fiber may transmit up to 1 liter/min of gas flow, cooling the fiber as well as clearing the surgical field of smoke or blood.

The system consists of following components:

- OmniGuide Surgical fiber-enabled laser system
- Medical Grade Helium gas connected to the laser system
- OmniGuide Surgical CO₂ Laser Fiber, sterile, single use
- OmniGuide Surgical handpiece or other rigid delivery device
- Gas Filter Unit (GFU) system for use of sterile gas in sterile field

SECTION II: INDICATIONS FOR USE

The OmniGuide BeamPath OTO-M Fiber is indicated for the incision, excision, ablation, vaporization and coagulation of body soft tissues.

SECTION III: WARNINGS

- Pressurized gas exiting the fiber tip during the laser procedure may cause gas embolism if fiber comes in contact with blood vessels or vascular tissue, leading to serious injury or death.
 - Do not bring the fiber tip into direct contact with blood vessels or vascular tissues
 - Monitor vital signs of the patient and blood oxygen level for symptoms of gas embolism
 - Stop using the laser if you see sudden or unexplained changes of the levels of end tidal carbon dioxide, subcutaneous emphysema, sudden change in oxygen tension via oximetry or sudden changes in hemodynamics
- Pressurized gas exiting the fiber's tip during the laser procedure may cause temporary inflation and separation of sub-mucosal flaps or emphysema under superficial layers of tissue.
 - To reduce this risk, use a longer working distance between the fiber tip and tissue.
 - Any signs of subcutaneous emphysema or hemodynamic instability during use, immediately stop using laser and proceed to standard treatment protocol.
- Inspect fiber packaging. Do not use if package is damaged or opened.
- Fiber is intended for single use only. Do not re-sterilize or reuse.
- Fiber reuse may pose the risk of infection.
- Re-sterilization/reprocessing of the fiber could lead to fiber damage, resulting in fiber failure leading to potential:
 - Laser radiation exposure
 - Fiber overheating
 - Fiber foreign material being left in the patient
- Handle the fiber with care. Do not bend it excessively (radius < 2 cm) to avoid damage and increased possibility of failure during use.
- Medical Grade Helium gas must be flowing through the fiber at specified settings throughout the procedure to avoid increased possibility of failure during use.
- The laser should be activated only when the laser treatment site is clearly visible and the fiber tip is directed at it.


- Laser energy can cause eye damage, burns, and ignite flammable materials.
 - User should wear appropriate laser eyewear and skin protection in accordance with the facility's laser safety policies
 - Avoid reflective materials and combustive materials near the surgical site
 - Flammable prep solutions must not be used during laser procedures
 - Fire extinguisher must be available for the operating room when the laser is being utilized
 - Refer to ANSI Z136.3 Use of Lasers in Healthcare Facility
- If working near open vestibule, lower pressure to 10psi and use laser in single pulse mode.

SECTION IV: PRECAUTIONS

-  Observe all laser and compressed gas safety precautions as specified in the laser system User Manual.
- Do not use the fiber outside its recommended parameters in Section VI. SPECIFICATIONS below. Always use the minimum appropriate laser power and gas settings required for a given procedure.
- See the Fiber Procedural Setting Guide (DOC-GEN-3) or the OR Staff Training Binder (DOC-FELS-27) for recommended power and gas ranges.
- Laser surgical procedures should be performed only by a licensed physician or dentist adequately trained and familiar with such surgical techniques and clinical use of CO₂ lasers.
- A smoke evacuation system or suction should be utilized to remove surgical plume.
- Do not touch the fiber tip to tissue or fluids as it will increase the risk of tip clogging and/or product failure. If tissue or debris buildup occurs, set gas pressure to at least 10psi and clean the fiber tip by dipping it in saline solution.
- There is a possibility that an OmniGuide BeamPath OTO-M Fiber will fail during the course of a surgical procedure. Therefore, do not start a surgical procedure unless you have at least two fibers or another surgical tool available.

SECTION V: ADVERSE REACTIONS

- Inflammatory reaction at the sites where energy is applied.
- As with any surgical tool, iatrogenic injury may occur.
- Pressurized gas may cause mechanical damage to the middle ear, inner ear anatomy and perilymph at pressures greater than 10psi.
 - Warning: reduce pressure to 10psi when working in the middle ear
- Continuous lasing without adequate cooling flow may cause the fiber tip to heat to temperatures > 60° C.
 - Warning: use fiber in single pulse mode in middle ear and near critical structures
 - Warning: do not use fiber in contact mode on critical structures and locations

SECTION VI: SPECIFICATIONS

		SPECIFICATIONS
LASER SETTINGS	Maximum input power (laser setting)	10 W (continuous or single pulse mode) Pause for 3 seconds between pulses
	Recommended minimum working distance between fiber tip and tissue	3 mm
	Medical grade helium input pressure	See the Fiber Procedural Setting Guide (DOC-GEN-3) or the OR Staff Training Binder (DOC-FELS-27) for recommended power and gas ranges.
	Gas flow at recommended pressure	<1 lpm

		SPECIFICATIONS
PRODUCT SPECS.	Fiber Length	1.5 m
	Nominal spot size at fiber tip	0.25 mm
	Beam divergence	200 mrad (11.46°, full angle)
	Maximum outer diameter of the distal portion of the fiber assembly	0.55 mm

COMPATIBLE LASER SYSTEMS AND SURGICAL INSTRUMENTS

OmniGuide fibers are only intended to be used with OmniGuide handpieces and lasers. Please see the enclosed fiber compatibility chart for OmniGuide fiber/handpiece compatibility. Please note that this chart may have changed between when you purchased from OmniGuide, and when you use your fiber or handpiece. If there is any discrepancy, the most up-to-date version of the fiber compatibility chart can be found at omniguideurgical.com/fibercompatibility.

DIRECTIONS FOR USE

 See laser system operating manual for detailed instructions, maintenance, and additional warnings and precautions.

 For sterile cases, use of a Gas Filter Unit (GFU) is required. See GFU instructions for details.

 Inspect fiber packaging. Do not use if package is damaged or opened.
Note: See product packaging for expiration date.

1. Use sterile technique to open the fiber package.

 Handle the fiber with care. Do not bend it excessively (radius <2 cm) to avoid damage and increased possibility of failure during use

2. Insert the fiber into a compatible surgical instrument

Note: See surgical instrument instructions for use for detailed instructions, additional warnings and precautions.

3. Attach the fiber connector to the laser adapter. Ensure fiber connector is securely locked to the adapter. Verify gas is flowing from the fiber tip.

4. Turn on gas supply and set gas pressure.

5. Activate the laser energy to cut, ablate, or coagulate soft tissue for the desired application.

 Do not use the fiber outside its recommended parameters in the Specification Table. Always use the minimum appropriate laser power and gas settings required for a given procedure.

 See the Fiber Procedural Setting guide (DOC-GEN-3) or the OR Staff Training Binder (DOC-FELS-27) for recommended power and gas settings

Note: The CO₂ laser beam diverges as it exits the distal end of the fiber. Refer to Specification Table for spot size at the fiber tip and beam divergence.

6. Turn off the laser system and gas supply after procedure completion.

7. Disconnect the fiber connector from the laser adapter and remove the fiber from the surgical instrument.

8. Visually inspect the fiber and verify that no portion is missing.

9. Treat and dispose of used fibers as biohazardous material according to facility guidelines.

FOR FURTHER INFORMATION

If further information on this product is needed, please contact OmniGuide Customer Service at (888) OMNI-GUIDE (888-666-4484) in the US, or your authorized distributor. See www.omni-guide.com/patents for information about the coverage of this OmniGuide® product by U.S. and/or foreign patent rights.

USA

OmniGuide, Inc.
4 Maguire Rd,
Lexington, MA 02421

Phone: +1 (617) 551-8444 | (888) OMNI-GUIDE

Fax: +1 (617) 551-8445

Email: sales@omni-guide.com | URL: www.omni-guide.com

EUROPEAN AUTHORIZED REPRESENTATIVE:

MediMark Europe, Sarl
11, rue Emile Zola - BP 2332
38033 Grenoble Cedex 2 France

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