
FMX™

- EN** FMX Reusable Power Module - Directions for Use
Must be sterilized prior to initial use
- FR** Module d'alimentation réutilisable FMX - Mode d'emploi
À stériliser avant la première utilisation
- DE** FMX-Wiederverwendbares Strommodul - Gebrauchsanleitung
Muss vor der erstmaligen Verwendung sterilisiert werden
- IT** Modulo di alimentazione riutilizzabile FMX - Istruzioni per l'uso
Sterilizzare prima dell'uso iniziale
- ES** Módulo de potencia reutilizable FMX - Instrucciones de uso
Debe esterilizarse antes de usarse por primera vez
- PT** Módulo de Alimentação Reutilizável FMX - Instruções de utilização
Tem de ser esterilizado antes da utilização inicial
- NL** Herbruikbare FMX-energiemodule - Instructies voor gebruik
Moet vóór het eerste gebruik worden gesteriliseerd

Device Description: The FMX Reusable Power Module (RPM) controls the delivery of power from the FMX Generator to an attached FMX surgical instrument. The RPM is intended for use exclusively with the FMX Generator and compatible FMX surgical instruments.


The RPM is a reusable device with a limited service life. The RPM has an internal counter that limits its service life to 100 procedures. The generator will notify the user after 100 procedures are completed. A completed procedure is registered in the RPM when it is unplugged from the generator and remains inactive for a period of time. The number of remaining procedures is indicated on the generator display when the RPM is connected to the generator.

Indications: The FMX Reusable Power Module (RPM) delivers the necessary power to operate various FMX surgical instruments. For specific indications of the individual FMX surgical instruments, see their respective Directions for Use.

Warranty: Please refer to the FMX Generator Directions for Use for warranty information.

Description of Symbols: The following symbols are found on the RPM labeling.

Rx ONLY
Prescription
Only

Warnings

The RPM requires cleaning and sterilization prior to initial use.

Devices labeled “for single use only” must not be reprocessed for re-use.

Single-patient use devices must not be reused, as they are not designed to perform as intended after the first usage. Changes in mechanical, physical or chemical characteristics introduced under conditions of repeated use, cleaning and re-sterilization may compromise the integrity of the device leading to diminished safety, performance and/or compliance with relevant specifications and standards.

The RPM should be thoroughly cleaned and sterilized between each use. FMX surgical instruments are single-patient use only and should never be reused.


Because the contamination of this device may include human blood or body fluids capable of transmitting bloodborne pathogens, the requirements of OSHA Standard 29 CFR 1910.1030 – *Occupational Exposure to Bloodborne Pathogens* should be followed.

Never use metal brushes or steel wool for cleaning the RPM.

Pre-Cleaning


- Detach the cooling air line from both ends of the RPM and discard (figure 1).
- Wipe the RPM with a sponge or a lint free towel moistened with water (not saline) to remove soil and visible contaminants.
- Keep the RPM moist by covering with a damp cloth until the cleaning process can take place.

The RPM requires cleaning before sterilization. **If manually cleaning the RPM**, follow the instructions below. **If using an automatic cleaning system**, skip to the next section.

Manual Cleaning Instructions


- Rinse the RPM with cool tap water and wipe all gross visible contamination with a soft towel.
- Completely submerge and soak the RPM in a pH neutral enzymatic detergent (prepared according to the manufacturer's instructions) for 55 - 60 seconds.
- Flush the air-line connection lumen with detergent solution using a 10 ml syringe a minimum of 3 times or until no visible soil or contaminants come out (figure 2).
- Manually clean the RPM with a soft-bristled brush or equivalent. Using a 0.419 inch diameter brush, scrub the lumened areas of the RPM for 55 - 60 seconds.
- Thoroughly rinse the RPM with purified water for 55 - 60 seconds.
- Visually inspect the RPM for cleanliness. If any visible contamination remains, repeat these steps.

NOTE: The validation of this cleaning process was completed with Steris® Prolystica 2X Concentrate Enzymatic Detergent.

Automatic Cleaning Instructions

- Flush the cooling air line connection lumen on both ends of the RPM with water using a 10 ml syringe a minimum of 3 times or until no visible soil or contaminants come out (figure 2).
- Remove the syringe prior to the automated cleaning cycle.
- Rinse any debris off of the RPM with cool water for a minimum of one minute or until all loose contaminants are removed (figure 3). The temperature of the water should not exceed 30° - 35° C.
- Scrub the device using a damp towel and a soft-bristled brush with a 0.419 inch diameter for a minimum of 1 minute. Pay particular attention to the lumened areas.

- Perform a visual inspection of the RPM to assure that all visible contamination has successfully been removed. If not, repeat the pre-cleaning and cleaning steps.
- Place the RPM in an automatic cleaning container.
- Place in an automated washer/disinfector with the following (default) cycle:

Phase	Settings	
Pre-Wash	Circulation time (Min and Max)	0:15 -15:00
	Water type	Cold tap water
	Temperature	< 16° C
Wash 1	Circulation time (Min and Max)	1:00- 5:00
	Injection Range	¼ - 2oz
	Water type	Hot tap water
Wash 2	Temperature	≥ 43° C
	Circulation time (Min and Max)	1:00 - 15:00
	Injection Range	¼ - 2oz
Rinse 1	Water type	Hot tap water
	Temperature	≥ 43° C
	Circulation time (Min and Max)	00:15 - 10:00
Thermal Rinse	Water type	Hot tap water
	Temperature	≥ 82° C
	Circulation time (Min and Max)	00:15 - 10:00
Pure Water Rinse	Water type	Cold tap water
	Temperature	< 16° C
	Circulation time (Min and Max)	00:15 (non-recirculation)
Drying	Circulation time (Min and Max)	00:00-60:00
	Water type	NO WATER
	Temperature	66° -115.5° C

- After the automated wash cycle is complete, blow any excess water out of the connectors with compressed air.

NOTE: The validation of this cleaning process was completed with Steris® Prolystica 2X Concentrate Enzymatic Detergent and a HAMO LS-1000 Washer-Disinfecter.

NOTE: Dilute detergent according to manufacturer's instructions.

After the drying cycle is complete, the product is ready for steam (autoclave) sterilization.

Steam Sterilization (Moist Heat) is Recommended




NOTE: An autoclave cycle has been validated as being capable of sterilizing this medical device; however autoclave design and performance can affect the efficacy of the process. Healthcare facilities should validate the process that they use, employing the actual equipment and operators that routinely process the devices.

NOTE: Autoclaves should comply with the requirements of, and be validated, maintained and checked in accordance with EN285/EN 13060, ISO 17665 and ANSI/AAMI ST79.

Preparing the RPM for Moist-Heat (Autoclave) Sterilization Process

Wrap the RPM with CSR wrap (or equivalent). Use only wrappers validated for use in double simultaneous wrapping. Two single-layer wrappers or one bonded double-layer wrapper can be used. Simultaneous double wrapping refers to when two layers of wrap material are wrapped together simultaneously using the fold technique described below.

- Place the wrap sheets on the table to form a diamond shape. Place the RPM to be wrapped in the center of the wrap, parallel with the edge of the table (figure 4).

- Bring the lower corner up to completely cover the contents and fold the tip back on itself to form a tab or flap (which is used later to assist in opening the pack aseptically) (figure 5).
- Fold the left corner over the contents and fold the tip back to form a tab (figure 6).
- Fold the right corner over the left fold and fold the tip back on itself to form a tab (figure 7).
- Bring the top corner down over the contents and tuck the corner under the right and left folds. A small tab should be left out for easy opening. Secure with steam indicator tape (figure 8).

Moist-Heat (Autoclave) Sterilization Process

- Pre-inspect and set up the autoclave according to manufacturer's instructions.
- Place the wrapped RPM into the steam sterilizer.
- Verify the following autoclave cycle parameters:

	US	EU
Method	Moist Heat Sterilization	Moist Heat Sterilization
Cycle	Pre-vacuum	Pre-vacuum
Temperature	270° F (132° C)	273.2° F (134° C)
Exposure Time	4 minutes	3 minutes
Drying Time	20-30 minutes	20-30 minutes

- After the sterilization cycle is complete, allow the RPM to cool before handling.
- Product is sterile and ready for use.

For detailed instructions on how to install and use the RPM with FMX surgical instruments, refer to the Directions for Use information packaged with each individual FMX surgical instrument.

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For customer support inquiries,
please contact us at 888.666.4484

Domain Surgical, Inc.
1370 South 2100 East
Salt Lake City, UT 84108
888.666.4484
omni-guide.com

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Domain Surgical and FMX are registered trademarks of Domain Surgical, Inc. All content is copyright of Domain Surgical, Inc. © 2017. Visit www.omni-guide.com/patents for information about the coverage of these products by U.S. and/or foreign patent rights. LB-037 Rev. E

 888.666.4484 | 617.551.8444
 customerservice@omni-guide.com
 www.omni-guide.com