
Instructions for Use - Single-use laser fibers

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R_X ONLY

Federal law restricts this device to sale by or on the order of a physician.

Content

1	About these instructions for use	5
1.1	Signal words used in these instructions for use	5
1.2	Precautions and safety information	5
2	Product information	7
2.1	Product description	7
2.2	Classifications	8
2.3	Patient Population	8
2.4	Contraindications	8
2.5	Potential side effects	8
2.6	Sterile on delivery	8
2.7	Laser fiber to laser device compatibility	9
2.8	Technical data of laser fibers	10
2.9	Compatible instrument / endoscope	11
2.10	Symbols used and their meaning	12
3	Use of the laser fiber	13
3.1	Transport and storage conditions	13
3.2	Checking the single-use laser fiber before use	14
3.3	Unpack and prepare the single-use laser fiber for surgery – step by step	14
3.4	Disposal	15

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1 About these instructions for use

These instructions for use provide important information concerning the safe operation of single-use laser fibers manufactured by LISA Laser Products GmbH, Germany (LISA Laser Products) as listed on the front cover.

Before using the laser fiber read these instructions for use and those for the laser device thoroughly! Follow the given instructions!

The safety instructions in this document are intended to prevent possible injuries to the patient, user or third parties and material damage.

Please keep these instructions for use for future reference.

1.1 Signal words used in these instructions for use

The following table shows the meaning and classification of the signal word to the left of the safety instruction.


Indicates a highly hazardous situation which if not avoided will result in death or serious injury.

Indicates a hazardous situation which if not avoided could result in death or serious injury.

Indicates a fairly hazardous situation which if not avoided could result in minor or moderate injury.

Indicates imminent material damage and suggests ways to avoid possible material damage.

Provides recommendations, information and advice for efficient use.

1.2 Precautions and safety information

The laser system (including laser device, laser fiber and compatible surgical instrument) may only be used by qualified personnel who have appropriate medical knowledge or professional clinical experience and have been trained in the correct and safe operation of the device based on the instruction for use.


Please also read and observe the precautions and safety information in the other chapters.

Any serious incident (i.e. the death of a patient, user or other person, the temporary or permanent serious deterioration of a patient's, user's or other person's state of health or a serious public health threat) occurring in relation to the device must be reported immediately to LISA Laser Products and the national competent authority responsible for medical device surveillance.

**Laser radiation**

Laser radiation can cause injury to the eyes and skin.

- Observe the necessary safety measures for laser application.
- Avoid direct laser radiation to the eyes and skin.
- Do not look into the laser beam.
- Avoid scattered radiation, which for example may arise from reflection on reflective surfaces of instruments.
- Wear laser safety goggles adequate for the radiation.

**Detachment of fragments of a damaged laser fibers.**

Fragments can detach from damaged laser fibers. Unretrieved device fragments (UDF) can cause injuries or diseases in the human body.


- Use only fiber fixations which exclude mechanical damage to the laser fiber.
- Use only instruments / endoscopes without sharp corners and edges. Check the instrument / endoscope before use.
- In combination with the laser fiber, use only compatible instruments / endoscopes. Compare the technical data such as bend radius and necessary working channel of the laser fiber and the provided instruments / endoscopes.

Rough handling of the laser fiber can lead to mechanical damage and fiber breakage.

Uncontrollable laser radiation can escape from the damaged laser fiber and cause burns and tissue damage.

- Handle the product carefully to avoid kinks and other damage to the laser fiber and / or to the fiber connector.
- Check the integrity of the laser fiber before use.

Rough handling during application or transportation of the laser fiber can lead to mechanical damage and fiber breakage.

Damaged laser fibers may not be used and may lead to a termination or delay of the procedure.

- Handle carefully to avoid kinks and other damage to the laser fiber or the fiber connector.
- Check the integrity of the laser fiber before use.
- Keep suitable replacement ready.

2 Product information

2.1 Product description

The described single-use laser fibers are products of LISA Laser Products and have been developed and tested for the use in combination with medical laser devices manufactured by LISA Laser Products.

See Table 1 for laser fiber – laser device - compatibility information and maximum laser power settings allowed.

The jacket of the laser fibers is made of high-quality plastics and protects the fiber against mechanical damage (Figure 1).


Figure 1: Laser fiber design

The terms used in Figure 1 for the elements of a laser fiber are not standardized and may differ from other publications.

The proximal end of the laser fiber material is bonded into a LISA F-SMA type fiber connector (Figure 2), with which the laser fiber is connected to the fiber port of the laser device.

The distal tip of the laser fibers described in these instructions for use are bare ended. The laser radiation emerging from the bare ended distal tip is diverging at a full cone angle of app. 26°.


Figure 2: LISA Laser Products single-use laser fiber

2.2 Classifications

The single use laser fibers are a 'medical device' according to the definition of Regulation (EU) 2017/745 Art. 2 (1) intended to be used for human beings for the treatment of disease.

Medical device class according to Council Directive 93/42/EEC (MDD)	IIb
Medical device class according to Regulation (EU) 2017/745 (MDR)	IIb
Medical device nomenclature according to EMDN	Z12011080
Medical device nomenclature according to GMDN	62600
Basic UDI-DI	42503419LaserFiberSU44

2.3 Patient Population:

Lisa laser fibers have not been tested and are not claimed to be included for use on vulnerable populations, including but not limited to children, pediatric, geriatric, pregnant/lactating women, immunocompromised individuals, elderly subjects, sepsis, anticoagulation therapy or bleeding disorders, or patients with implants, and risk of reproductive effects remain unknown. However, the state-of-the-art, has demonstrated the use of Laser and fibers to be safe and effective when utilized by a qualified physician in a qualified clinical setting.

2.4 Contraindications

The laser fibres are not intended for direct application on the central nervous system and on the cardiovascular system.

Contraindications that are directly linked with the laser fibres are not known. The contraindications of surgical and endoscopic laser applications are the same as those of conventional endoscopy or surgery in the concerned field.

Contraindications to receive therapeutic laser procedures are a general inability to receive surgical or endoscopic procedures, pregnancy, sepsis, anticoagulation therapy or bleeding disorders. In each case the user must decide on the basis of the patient's condition whether the application of the laser system is appropriate

2.5 Intended use

Transmission of laser radiation during therapeutical laser procedures. The laser fibers are intended for transient use.

2.6 Contraindications

The fibers are not intended for direct application on the central nervous system and on the cardiovascular system.

Contraindications that are directly linked with the laser fibers are not known. The contraindications of surgical and endoscopic laser applications are generally the same as those of conventional endoscopy or surgery in the concerned field. Contraindications are a general inability to receive surgical or endoscopic procedures, pregnancy, sepsis, anticoagulation therapy or bleeding disorders. In each case the user must decide on the basis of the patient's condition whether the application of the laser system is appropriate.

2.7 Potential side effects

Potential complications which may result from a surgical and endoscopic laser procedure are similar to those normally encountered in conventional surgery in the same field and may include but are not limited to: thermal damage to surrounding structures, pain, perforation, haemorrhage, infections, sepsis, tissue trauma, bleeding and oedema.

2.8 Sterile on delivery

The single-use laser fibers are delivered EO sterilized in separate see-through pouches and delivered in packs of five per pack.

The laser fibers are not made for reprocessing and reuse.

Re-processing of single-use laser fibers may cause damage to the laser fiber.

Fragments can detach from damaged laser fibers. Unretrieved device fragments (UDF) can cause injuries or diseases in the human body.

- Do not re-process the laser fiber.
- Use the laser fiber not more than once.

Re-processing of single-use laser fibers

The reprocessing of laser fibers for single use may result in infections by insufficient sterility.

- Do not re-sterilize the laser fiber.
- Use the laser fiber not more than once.

2.9 Laser fiber to laser device compatibility

LISA Laser Products laser fibers differ in optical properties, in design, in optical core and outer diameter. Table 1 shows the compatibility of LISA Laser Products single-use laser fibers with LISA Laser Products laser devices and the maximum permissible power settings.

Table 1: Laser fiber to laser device compatibility

Laser fiber Order number	 Laser device compatibility¹						
	<i>Sphinx</i>	<i>Sphinx litho</i>	<i>Sphinx jr.</i>	<i>RevoLix 200</i> <i>RevoLix 120</i>	<i>RevoLix DUO</i>	<i>RevoLix jr.</i>	<i>RevoLix HTL</i> <i>RevoLix HTL cw</i> <i>RevoLix HTL eco</i>
<i>LithoFib-SU</i> 101 503 576	-	✓ 12 W and max. 15 Hz	✓ ⁴ 12 W and max. 2.0 J and max. 14 Hz	-	-	✓	✓ 20 W
<i>FlexiFib-SU</i> 101 503 387	✓ ² 30 W and max. 3.5 J	✓ ² 30 W and max. 3.5 J	✓ ² 30 W and max. 3.5 J	✓ 40 W	✓ ³ 40 W 20 W	✓	✓ 40 W
<i>SureFib-SU</i> 101 503 513 101 503 648	-	✓ ² 30 W and max. 4.0 J	✓ ² 30 W and max. 3.5 J	✓ ⁵ 40 W	-	✓	✓ 40 W
<i>SureFib-SU 5m</i> 101 503 569	-	✓ ² 30 W and max. 4.0 J	✓ ² 30 W and max. 3.5 J	✓ ⁵ 40 W	-	✓	✓ 40 W
<i>PercuFib-SU</i> 101 503 384	✓ ² 50 W and max. 3.5 J	✓	✓	✓ 50 W	✓ ³ 50 W 20 W	✓	✓ 50 W
<i>RigiFib-SU</i> 101 503 289 101 503 645	✓	✓	✓	✓	✓	✓	✓ 150 W
<i>RigiFib-SU 800</i> 101 503 385	✓	✓	✓	✓	✓	✓	✓ 150 W
<i>RigiFib-SU 1000</i> 101 503 386	✓	✓	✓	✓	✓	✓	✓ 150 W


¹ - not compatible with laser, ✓ compatible with laser

² maximum power and maximum energy (meet both limits)

³ upper value: max. power cw laser (Tissue Mode),
lower value: max. power pulsed laser (Stone Mode)

⁴ providing *LithoFib* compatibility is shown on the start display and *LithoFib* is selected

⁵ Not for use with RevoLix 50, RevoLix 70 or power upgraded devices

**CAUTION**

Absorption of laser power in the fiber connector by combining the laser fiber with incompatible laser devices

Absorption of laser power in the fiber connector may heat up the fiber connector and may cause burns when touching the hot fiber connector.

- Comply with the specified compatibility requirements.

NOTICE

Mechanical damage of the laser fiber caused by combining it with incompatible instruments / endoscopes or laser devices.

Using incompatible instruments / endoscopes or laser devices can result in loss of function of the laser fiber.

- Use only compatible instruments / endoscopes in combination with the laser fiber. The bend radius / curvature of the instrument / endoscope must not be less than the minimal bend radius of the laser fiber. The outer diameter of the laser fiber is required to be smaller than the working channel diameter of the instrument / endoscope.
- Comply with the specified compatibility requirements.
- Use only instruments / endoscopes without sharp corners, edges and / or kinked working channels. Check the instrument / endoscope before use.
- Use only fiber fixations which do not damage the laser fiber.
- Insert the laser fiber carefully into the instrument / endoscope.

2.10 Technical data of laser fibers

Technical data of LISA Laser Products single-use laser fibers are listed in Table 2.

Table 2: LISA Laser Products single-use laser fibers technical data

Laser fiber Order number	Optical core diameter	Outer diameter	Length	Min. bend radius during operation
<i>LithoFib-SU</i> 101 503 576	200 µm	500 µm / 1.5 Fr	2.5 m	≥ 10 mm
<i>FlexiFib-SU</i> 101 503 387	272 µm	420 µm / 1.3 Fr	2.5 m	≥ 15 mm
<i>SureFib-SU</i> 101 503 513 101 503 648	272 µm	420 µm / 1.3 Fr	2.5 m	≥ 10 mm
<i>SureFib-SU 5m</i> 101 503 569	272 µm	420 µm / 1.3 Fr	5.0 m	≥ 10 mm
<i>PercuFib-SU</i> 101 503 384	365 µm	730 µm / 2.2 Fr	2.5 m	≥ 40 mm
<i>RigiFib-SU</i> 101 503 289 101 503 645	550 µm	750 µm / 2.3 Fr	2.5 m	≥ 60 mm
<i>RigiFib-SU 800</i> 101 503 385	800 µm	1200 µm / 3.6 Fr	2.5 m	≥ 70 mm
<i>RigiFib-SU 1000</i> 101 503 386	940 µm	1400 µm / 4.2 Fr	2.5 m	≥ 120 mm

2.11 Compatible instrument / endoscope

Any instrument / endoscope is deemed to be compatible with a laser fiber providing the instrument's / endoscope's working channel is

- curved at a bend radius equal or larger than tolerated by the laser fiber to be used,
- free from sharp corners and edges (including the fiber fixation),
- larger in diameter than the outer diameter of the laser fiber to be used.

A compatible instrument / endoscope will accept the insertion of a laser fiber without any force.

Emission of laser radiation caused by mechanical damage of the laser fiber. Laser fibers can be damaged by combination with incompatible instruments / endoscopes or laser devices.

Fiber breakage may cause unintended emission of laser radiation from the damaged laser fiber and may result in heat, tissue damage and burns.

- **Use the laser fiber only in combination with compatible instruments / endoscopes. Compare technical data of the laser fiber and the instruments / endoscopes provided, especially working channel and bend radius, prior to use.**
- **Comply with the specified compatibility requirements.**
- **Use only instruments / endoscopes without rough edges. Check the instrument / endoscope prior to use.**
- **Use only laser fiber fixations which prevent mechanical damage to the laser fiber.**
- **Insert the laser fiber carefully into the instrument / endoscope.**

Insufficient fixing of the laser fiber.

If insufficiently fixed the laser fiber can slip out of position into the instrument / endoscope and out of sight of the surgeon. Laser radiation can cause heating of the instrument / endoscope and tissue damage and burns to the patient.

- **Fix the laser fiber on / in the instrument / endoscope.**
- **Only use fixings that are suitable for the outside diameter of the laser fiber. Compare the specifications of the laser fiber (outer diameter) and the intended fixation.**

Excessive bending of a laser fiber.

If bent at a smaller bend radius than specified, laser radiation can escape at the bend. The resulting heat-development can cause tissue damage and burns.

- **Do not fall below the minimum bend radius.**

2.12 Symbols used and their meaning

Quantity	Batch code	Order number	Do not re-use	Do not re-sterilize	Medical device
Do not use if package is damaged	Single sterile barrier system with protective packaging inside	Single sterile barrier system	Outer diameter (Jacket)	Core diameter (Optical core)	Length
Date of manufacture (year-month)	Use-by date (year-month)	Be aware of warnings and safety precautions	Consult instruction for use	Latex free	Sterilized using ethylene oxide
		Made in Germany			
Minimum bend radius	Manufacturer	Certificate of origin	CE-Mark	US federal law restricts this device to sale by or on the order of a physician!	
Temperature limitation	Keep away from sunlight	Keep dry	Humidity limitation	Atmospheric pressure limitation	

The product label (attached to the sterile barrier system and to the package) comprises all necessary data for the identification of the laser fiber as well as safety relevant information.


Figure 3: Product label of a **SureFib-SU** (example)

3 Use of the laser fiber

3.1 Transport and storage conditions

Table 3: Transport and storage conditions

	Temperature	Relative humidity	Atmospheric pressure
Transport	-18 °C to +60 °C 0 °F to 140 °F	30 % to 85 %	700 hPa to 1060 hPa
Storage	+5 °C to +30 °C 40 °F to 86 °F	30 % to 75 %	700 hPa to 1060 hPa

Store the laser fiber in a dry place protected against direct light, heat radiation and dust.


WARNING

Rough handling of the sterile package, inadequate transport and storage conditions can lead to mechanical damage of the sterile barrier of the package.

A damaged sterile package can cause infection by insufficient sterility.

- Check the sterile package before each use, transport or storage for the integrity of the sterile barrier.
- Do not use a product whose sterile package is damaged. Dispose off the product immediately in case of suspected damage to avoid the risk of infection.
- Do not use sharp or pointed objects in handling the sterile package.
- Observe the specified transport and storage conditions.


CAUTION

Mechanical and climatological stress during transport / storage of the laser fiber.

Mechanical damage and fiber breakage.

- Comply with the transport and storage conditions.
- Check the integrity of the laser fiber and package before use.

CAUTION

Rough handling in the application, processing or transportation of the laser fiber can lead to mechanical damage and fiber breakage.

Sharp edges can cause injury.

- Handle carefully to avoid kinks in the laser fiber and other damage to the laser fiber or the fiber connector.
- Check the integrity of the laser fiber before use.


3.2 Checking the single-use laser fiber before use

1. Check the product label for legibility before using the laser fiber. Do not use the laser fiber, if the information is unreadable or unavailable.
2. Inspect the pouch for any breach or puncture of the package to ensure sterility. If the package is damaged or sterility is in question, do not use the product.
3. Make sure that the use-by date is not exceeded.
4. Inspect the laser fiber for kinks, fractures or other damage. The laser fiber has to be free of damage and kinks on its entire length.

3.3 Unpack and prepare the single-use laser fiber for surgery – step by step

1. The o.r. scrub nurse removes the sterile laser fiber from the pouch and inspects the laser fiber for kinks, fractures or other damage. Test the mechanical integrity of the distal tip by applying a sideways force as doing a full stop with a ballpoint pen (Figure 4).

Figure 4: Test of mechanical integrity

2. The o.r. scrub nurse hands over the fiber connector to the circulating nurse, who operates the laser system.
3. The circulating nurse removes the protective cap from the fiber connector.
4. Connect the laser fiber to the laser device. Refer to instructions for use of the laser device for detailed information how to connect the laser fiber.
5. Set the aiming beam to a high intensity. Check if the aiming beam is emitted from the distal tip of the laser fiber and not from the length of the laser fiber. A proper laser fiber projects a circular image of the aiming beam on a white non- reflective target.
6. If the aiming beam is weak or poorly visible or appears somewhere along the length of the laser fiber or from the stripped length at the distal tip, do not use the laser fiber.
7. Insert the distal tip of laser fiber into the compatible instrument / endoscope (see 2.9).
8. The distal tip of the laser fiber burns back slightly during surgery. The burn back may alter the optical beam quality of the laser radiation emerging from the distal tip and may change the surgical effect. In any case the laser fiber tip has to protrude at least 1 - 3 mm from the instrument / endoscope.

Emission of laser radiation with the distal fiber tip within the instrument / endoscope.

Emitted laser radiation may cause tissue damage and burns.

- **Never activate the laser unless the distal fiber tip protrudes from the instrument / endoscope by 1 mm at least.**

3.4 Disposal

The laser fiber gets in contact with tissue and body fluids and after surgery may present a potential source for infection.

Dispose the single-use laser fiber immediately after use in an unbreakable semi-solid container which encloses the disposed laser fiber safely. Follow the regulatory, technical and organizational required safety precautions for the disposal of septic, pointed and breakable items.

After use laser fibers are contaminated with biological material.

Contact with a used laser fiber may result in infections.

- **For the disposal of the laser fiber use containers that are suitable for biologically contaminated material or clean and disinfect the laser fiber before discarding it.**
- **Observe national regulations for disposal.**

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