



NexScope

Nephroscope and Mini Nephroscope




Key Features

- Outstanding performance: contrast enhanced image bundle, best vision
- Safe and reliable: robust, tightly sealed, laser welded, easy to use and to clean
- Direction of view ideally aligned at 8°
- Strong irrigation flow enables high stone-free rates
- Long lifetime for lowest cost of ownership
- Lightweight design
- Detachable instrument bridge for optimized cleaning and inspection
- Compatible with most popular laser and lithotripter systems
- Autoclavable: 5 min. @ 134 °C
- Made in Germany



Technical Specifications **NexScope** Nephroscope and Mini Nephrosopes



NexScope

Description	Part number	
NexScope 19 F 220	working channel 10 F, field of view 84°	101 503 437
NexScope Sheath 24 F	with obturator, colour code: red	101 503 440
NexScope Cleaning Brush, Working Channel 10 F	colour code: blue	101 503 486

Mini NexScope

Mini NexScope 13.5 F 220, 220 mm length	working channel 6 F, field of view 84°	101 503 435
Mini NexScope 13.5 F 320, 320 mm length	working channel 6 F, field of view 84°	101 503 436
NexScope Sheath 18 F	with obturator, colour code: yellow	101 503 439
NexScope Cleaning Brush, Working Channel 6 F	colour code: white	101 503 487

Consumables and Accessories for NexScope and Mini NexScope

NexScope Instrument bridge	for NexScope and Mini NexScope	101 503 438
Cleaning Brush for Side Port*	of instrument bridge, colour code: white	101 503 488
Cleaning Brush for Main Port*	of instrument bridge, colour code: blue	101 503 489
Three Slot Membrane	disposable, 1 package à 10 pieces	101 503 630
Sealing caps, for fibre fixation	disposable, 1 package à 10 pieces	101 503 499
Stop cock plugs*	1 package à 10 pieces	101 503 509
Lubricant PARALIQ GB 363*	40 g	101 503 634

* optional



Specifications are subject to change without notice.
Made in Germany 2019 - 02
Brochure Order no. 036 013 020

CE acc. Council Directive 93/42/EEC

U.S. federal law restricts these devices to sale by or on the order of a physician



LISA Laser Products GmbH
Albert-Einstein-Str. 4
37191 Katlenburg-Lindau
Germany
fon: +49 5556 9938-0
fax: +49 5556 9938-10
info@lisalaser.de www.lisalaser.de

